

DPST Accessories and Auxiliary

Our help- and needful DPST accessories and auxiliary is carefully adjusted to the demand of the DPST electrodes and also qualified for general laboratory praxis.

DPST Reference and Salt Bridge Solutions

DPST reference and salt bridge solutions are provided in 100 ml plastic containers. They are useful for the maintenance of DPST reference electrodes, the adjustment of the concentration of the reference electrolyte and for salt bridges.

DPST Reference and Salt Bridge Solution	Compounds
Electrolyte	KCl, c = sat. (saturated)
	KCl, c = sat., AgCl, c = sat.
	KCl, c = 3,0 mol/l
	KCl, c = 3,0 mol/l, AgCl, c = sat.
	KCl, c = 1,0 mol/l
	KCl, c = 1,0 mol/l, AgCl, c = sat.
	KCl, c = 0,1 mol/l
	KCl, c = 0,1 mol/l, AgCl, c = sat.
	K ₂ SO ₄ , c = sat.
	KNO ₃ , c = sat.
	KNO ₃ , c = 3,0 mol/l
	KNO ₃ , c = 1,0 mol/l
	KNO ₃ , c = 0,1 mol/l
	(NH ₄) ₂ SO ₄ , c = 1,0 mol/l
	(NH ₄) ₂ SO ₄ , c = 0,1 mol/l
	NaCl, c = 1,0 mol/l
	NaCl, c = 0,1 mol/l
LiOAc, c = 1,0 mol/l (lithium acetate)	

the company of innovations

DPST BEHNERT GMBH
DIREKT POTENTIOMETRISCHE SENSOR TECHNIK